



LASAK BioniQ QR Exocad library structure



Scanbody QR, indexed – long
Ref. No. 2815.00



bioniq_qr_premill_ag

bioniq_qr_premill_nt

Dummy base shown
in the Exocad library





**= CAM
replaced**

**= CAM
replaced**

BioniQ QR Premill abutment AG

BioniQ QR Premill abutment NT



The real base of the
BioniQ component



Ref. No.
2196.00

Ref. No.
2184.00

Ref. No.
2858.00



S3.5
S4.0
T4.0
S5.0
T5.0



exocad

LASAK BioniQ QN Exocad library structure



Scanbody QN, indexed – long
Ref. No. 2840.00



bioniq_qn_premill_ag

bioniq_qn_premill_nt

Dummy base shown
in the Exocad library





CAM
replaced



CAM
replaced

BioniQ QN Premill abutment AG


BioniQ QN Premill abutment NT




Ref. No.
2197.00

Ref. No.
2187.00

The real base of the
BioniQ component



Ref. No.
2859.00



S2.9